Gaumard®

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MONITOR SCENARIO CPR

Gaumard® Vitals will be enabled once connection is established.

Scenario Time

00:00:00

NIBP
Manual

Temperature

SpO2

MAP

RR

HR

LOG

00:03:38 Default Clinical State: Cardiac Rhythm Sinus;

00:00:34 On Demand: Cardiac Rhythm LBBB;

Team Add to log

sec

Apply
00:00

01:10:34

Quick Start HAL®

V

U

H

H

Power

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Gaumard® Vitals will be enabled once connection is established.

LOG

00:03:38 Default Clinical State: Cardiac Rhythm Sinus;

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ALL VITALS

Neurological 2

- Intracranial Pressure7.5 mmHg
Autoregulation: 100 %
- SeizuresNone

Cephalic 1

- Cyanosis0 %

Ocular 4

- Blink State15 blinks/min
- Pupil DilationRight: 5
Left: 5
- Pupil ReactionRight: On
Left: On
- Pupil Reaction Time0.5 sec

Airway 4

- Tongue EdemaOff
- LaryngospasmOff
- Pharyngeal SwellingOff
- Throat SoundNormal (3)

Breathing 14

sec

Apply00:00

MONITORSCENARIOCPR

ECG II

ABP

Heart Rate

ABP

MAP

Gaumard® Vitals will be enabled once connection is established.

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
Clinical State	Details	Adjust
1st deg AVB 85 bpm, BP 100/60		
Bradycardia RR 10, HR 45, 2nd deg. AVB (Type II)		
Christian supraventricular tachycardia		
Closed airway all difficult airway actuators on		
Evan sinus tach. with tachypnea		
Frank atrial fib. 180bpm. RR 22/min. BP 90/50. O2-Sat 93%		

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
The screenshot shows the 'Speech' window with several annotations in yellow boxes:

- 搜尋欄** (Search Bar): Points to the 'Search in Complaint' input field.
- 選單** (Menu): Points to the 'Libraries' list on the left, which includes 'Number', 'Occurrence', 'Reply', 'Urgency', and 'Custom'.
- 音量調節** (Volume Control): Points to the 'Speech Volume' slider at the bottom left, set to 15.
- 語音語言** (Speech Language): Points to the 'English' dropdown menu.
- 設定** (Settings): Points to the 'Settings' button in the top right.
- 語音重複時間** (Speech Repeat Time): Points to the 'Loop' checkbox and '0 sec' timer at the bottom right.

Libraries	Play Speech	Test
Number	Ankle broken	[Play]
Occurrence	Arm broken	[Play]
Reply	Blood in toilet	[Play]
Urgency	Can't catch breath	[Play]
	Can't move	[Play]
Custom	Can't move legs	[Play]

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




A form for creating a new group. It has a 'Name' field, a search bar, and a '+ New Provider' button. Below this, there are two sections: 'All Providers' and 'Group's Providers'. The 'Group's Providers' section has a green arrow button.

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

Search Speech

Language English

History	Asthma
Complaint	Diabetes
Number	Emphysema
Confusion	Had heart attack
Evasion	High blood pressure

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Settings

Details

Action Groups	Actions
Basic	Determine patency
Airway	Airway management
Breathing	Open airway
Circulation	Intubation Check
Trauma Care	Intubation evaluation
意識確認	Ventilate
	Extubate

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

A screenshot of a 'New Evaluation' form. The form has a dark background. At the top, it says 'New Evaluation'. Below this are two input fields: 'Name' and 'Description'. Both fields have red placeholder text. Below the 'Description' field, there are two sections: 'Categories' and 'Single entries'. Each section has a plus icon, a list of items (represented by up and down arrows), and a trash icon. At the bottom right of the form, there are two buttons: 'Cancel' and 'Save'.

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Adenosine

Concentrations: 2

Routes: Injection (INJECTABLE), Intravenous (SOLUTION)

Adrenalin

Concentrations: 1

Routes: Intravenous (SOLUTION), Intramuscular (SOLUTION), Subcutaneous (SOLUTION)

Albumin 5%

Concentrations: 0

Albuterol Sulfate

Concentrations: 1

Routes: Inhalation (SOLUTION), Oral (TABLET)

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

Concentrations / Syringes Routes

+ New Concentration / Syringe

Concentration	Description	Syringe	Amount	Edit	Delete
0.5 MG / 1 mL	0.5 MG / 1 mL		0.5		

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

A screenshot of a user profile form. On the left is a circular profile picture of a woman with the text 'Click on image to update it' below it. To the right of the image are input fields for 'Name' (containing '黃義哲'), 'Password' (masked with dots), 'Confirm Password' (masked with dots), 'Email' (containing 'fmost.ps@gmail.com'), and 'Phone' (containing '88622448545').

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The screenshot shows the 'Action' tab of a medical device interface. It features a search bar, a 'Selected Provider' dropdown, and a 'Settings' button. A list of 'Action Groups' (Basic, Airway, Breathing, Circulation) is on the left, and a table of 'Actions' is on the right. A yellow box highlights the 'Call for assistance' action, which has a green checkmark and a plus icon. A yellow box also highlights the 'Electrotherapy' section at the bottom, which includes 'SHOCK' and 'PACER' settings. Annotations in Chinese point to the 'Operator' (操作者) field, the 'Action Groups' (種類) list, and the 'Electrotherapy' (電擊) section.

Provider Action Evaluation

Actions 2

Search

Selected Provider Team Settings

Operator

Action Groups

- Basic
- Airway
- Breathing
- Circulation

Actions	Apply	Details
Assess responsiveness	✓	+
Call for assistance	✓	+

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
LOG	
04:40:00	[黃奕哲] Action: Establish IV
04:39:48	[奕哲] Action: Call for assistance
04:39:33	Action: Call for assistance
04:18:26	Action: Assess responsiveness

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Search


Filtered out: 0

Scenario	Adjust
Anthony 77-yr old man is found pulseless and apneic.	
Bradycardia Note: for this scenario to function as intended the instructor should enable automatic pacing capture in the 'Setup -> Auto-Responses'...	
Brent 55-year old male complains of substernal chest pain radiating to the jaw. Patient is diaphoretic, trembling, and has a waxing and waning...	
Christian 18-yr old college freshman calls EMS after taking a caffeine preparation prior to exams. He is short of breath and says that his heart is "b..."	
Darrell 65 year-old complains of dizziness and shortness of breath. He has a history of emphysema and does not know how much "breathing m..."	
Evan 72 year-old is found tachycardic and tachypneic. He is disoriented and barely responsive. Patient's wife says that he has both heart and l...	
Frank 50 year-old found at home complains of weakness and fatigue. He is conscious and alert but slightly short of breath. An empty bottle o...	
Gerard 35 year-old man calls for help at local restaurant after onset of alarming respiratory symptoms. (Patient experiencing allergic reaction to...	
Pulseless Arrest A young male was found unconscious.	
Tachycardia - pulse SVT with pulses. Note: this scenario can progress directly to the Pulseless Arrest scenario. To do so, load "Pulseless Arrest" after stopping...	

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Anthony

Back

Description : 77-yr old man is found pulseless and apneic.

Checklist

The checklist is empty

Items

- Anthony asystole 00:00:00 ⓘ
- Wait Indefinitely drug therapy expected --:--:-- ⓘ
- V Fib, coarse 00:00:00 ⓘ
- Wait Indefinitely electrical therapy expected - 300J threshold --:--:-- ⓘ
- V Tach, stable RR 16/min, HR 120bpm, BP 100/50 00:00:10 ⓘ
- Wait Indefinitely Hold in stable V. Tach --:--:-- ⓘ
- Healthy Resting normal vitals 00:00:20 ⓘ

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