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PALETTE

99% O2-Sat

Name	Description
1st deg AVB	85 bpm, BP 100/60
2nd Degree Heart Blo...	HR = 40, BP = 100/40, RR = 6
2nd Degree Heart Block	HR = 40, BP = 100/40, RR = 6
99% O2-Sat	
Anthony	asystole
apnea	RR 0/min
Asystole rhythm	patient is in asystole
Asystole	no electrical activity, apneic
Atropine	HR +10%, BP +7%
BP +40%	40% increase in blood pressure
BP -10%	10% reduction in blood pressure
BP -30%	30% reduction in blood pressure

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




Ventilation Duration:

Release Threshold: mmHg

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
COMPRESSION PROVIDER	VENTILATION PROVIDER	REPORTS
<input type="text" value="Henry"/>	<input type="text" value="Carl"/>	<input type="button" value="REPORTS"/>

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

The screenshot shows the LAB interface with two main sections: Lab Report and Lab Template. The Lab Report section has a table with columns: Name, Template, Description, and Sent. To the right of this table is a green button labeled 'New Lab Report' and a list of actions: View, Sort by Color, Selected Item, Send to monitor, Delete, Preview, and Edit... The Lab Template section has a table with columns: Template and Description. To the right of this table is a green button labeled 'New Template' and a list of actions: Selected Item, Edit..., and Delete. At the bottom of the interface is a navigation bar with tabs: QUICK LAUNCH, VIRTUAL MONITOR, FETAL MONITOR, PALETTE, SCENARIO, LABOR, LAB (selected), SPEECH, CPR, and PROVIDER ACTIONS.

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This close-up shows the 'New Lab Report' button in the LAB interface, which is highlighted with a red rectangular box. Below the button are the 'View', 'Sort by Color', and 'Selected Item' options.

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The screenshot shows the 'Edit Lab Report' window. It has a 'Save Options' section with four radio buttons: Healthy (selected), Critical, Care Required, and Other. Below these are 'Clear', 'Save', and 'Cancel' buttons. The 'Lab Report Information' section contains a 'Lab Template' dropdown menu with options: Arterial Blood Gas, CBC Differential (female), CBC Differential (male), and Glasgow Coma Scale. Other fields include 'Template Description', 'Lab Name', 'Description', 'Patient Name', 'Patient Age', 'Patient Gender', 'Time' (set to 00:00:00), and 'Report Notes'. At the bottom is a section for 'Lab Report Results'.

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Lab Report Information			
Lab Template:	<div><input checked="" type="checkbox"/> Arterial Blood Gas</div> <div><input checked="" type="checkbox"/> CBC Differential (female)</div> <div><input type="checkbox"/> CBC Differential (male)</div> <div><input type="checkbox"/> Glasgow Coma Scale</div>	Patient Name:	
Template Description:	Arterial blood gas analysis results, Complete Blood Count w/ Differen...	Patient Age:	
Lab Name:		Patient Gender:	
Description:		Time:	00:00:00 ▾
		Report Notes:	

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Edit Lab Report				
Save Options				
<input checked="" type="radio"/> Healthy	<input type="radio"/> Critical	<input type="radio"/> Care Required	<input type="radio"/> Other	
Clear	Save	Cancel		
Lab Report Information				
Lab Template:	<div><input checked="" type="checkbox"/> Arterial Blood Gas</div> <div><input checked="" type="checkbox"/> CBC Differential (female)</div> <div><input type="checkbox"/> CBC Differential (male)</div> <div><input type="checkbox"/> Glasgow Coma Scale</div>	Patient Name:	Noelle	
Template Description:	Arterial blood gas analysis results, Complete Blood Count w/ Differen...	Patient Age:	27	
Lab Name:	Report 1	Patient Gender:	female	
Description:	test	Time:	12:00:00 ▾	
		Report Notes:		
Lab Report Results				
Variable	Min Reference Value	Max Reference Value	Value	Comment
RBC	3.8	5.2	7	
Hgb	12	15	15	
Hct	35	45	40	
MCV	80	100	90	
MCH	26	35	28	
MCHC	32	36	35	
RDW	11.5	14.5	12	

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Lab Report			
Name	Template	Description	Sent
Report 1	Arterial Blood Gas, C...	test	

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Edit... Preview

Send to Virtual Monitor

Arterial Blood Gas, CBC Differential (female)

Report Time: 12:00:00
 Patient Name: Noelle
 Age: 27
 Sex: female
 Report Notes:

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SPEECH

Type

- ☐ PRERECORDED
- ☐ COMPLAINT
- ☐ CONFUSION
- ☐ DESCRIPTION
- ☐ EFFECT
- ☐ EVASION
- ☐ HISTORY
- ☐ LOCATION
- ☐ NUMBER
- ☐ OCCURRENCE
- ☐ REPLY
- ☐ URGENCY

Add Type

Options

Loop

☐ Loop 0 sec

Voice Recorder

REC

Output Gain

Voice Activation

Off 10

Language

☒ English ☐ Italian

☐ Finnish ☐ Spanish

☐ Russian ☐ French

QUICK LAUNCH
VIRTUAL MONITOR
PALETTE
SCENARIO
CPR
LAB
SPEECH
PROVIDER ACTIONS

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Type	COMPLAINT					
<input type="radio"/> PRERECORDED	<input type="radio"/> Ankle broken	<input type="radio"/> Arm broken	<input type="radio"/> Blood in toilet	<input type="radio"/> Can't catch breath	<input type="radio"/> Can't move	<input type="radio"/> Can't move legs
<input checked="" type="radio"/> COMPLAINT	<input type="radio"/> Chest hurts	<input type="radio"/> Coughing up blood	<input type="radio"/> Elephant on chest	<input type="radio"/> Feel dizzy	<input type="radio"/> Feel nauseous	<input type="radio"/> Feel weak
<input type="radio"/> CONFUSION	<input type="radio"/> Heart beating fast	<input type="radio"/> Heart pounding	<input type="radio"/> Heart trying to jump	<input type="radio"/> Hurt all over	<input type="radio"/> Hurts when breathing	<input type="radio"/> I've been cut
<input type="radio"/> DESCRIPTION	<input type="radio"/> Jaw hurts	<input type="radio"/> Left arm hurts	<input type="radio"/> Leg is broken	<input type="radio"/> Passing blood	<input type="radio"/> Peeing blood	<input type="radio"/> Pooping blood
<input type="radio"/> EFFECT	<input type="radio"/> Puking blood	<input type="radio"/> Short of breath	<input type="radio"/> Shoulder hurts	<input type="radio"/> Somebody shot me	<input type="radio"/> Stomach hurts	<input type="radio"/> Worst headache
<input type="radio"/> EVASION						
<input type="radio"/> HISTORY						
<input type="radio"/> LOCATION						

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